

NDEL1 Polyclonal Antibody

Cat No: HR1AP10115

For research use only

Overview

Product Name	NDEL1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	37kDa
GeneID?Human?	NDEL1 EOPA MITAP1 NUDEL
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	nudE neurodevelopment protein 1 like 1(NDEL1) Homo sapiens This gene encodes a coiled-coil protein that plays a role in multiple processes including cytoskeletal organization, cell signaling and neuron migration, outgrowth and maintenance. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene, and a pseudogene of this gene is located on the long arm of chromosome X. [provided by RefSeq, Mar 2012],