

UBP28 Polyclonal Antibody

Cat No: HR1AP10098

For research use only

Overview

Product Name	UBP28 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	118kDa
GeneID?Human?	USP28 KIAA1515
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	ubiquitin specific peptidase 28(USP28) Homo sapiens The protein encoded by this gene is a deubiquitinase involved in the DNA damage pathway and DNA damage-induced apoptosis. Overexpression of this gene is seen in several cancers. [provided by RefSeq, Oct 2016],