

# NTAL Polyclonal Antibody

Cat No: HR1AP10088

For research use only

## Overview

Product Name	NTAL Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	26kDa
GeneID?Human?	LAT2 LAB NTAL WBS15 WBSCR15 WBSCR5 HSPC046
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	linker for activation of T-cells family member 2(LAT2) Homo sapiens This gene is one of the contiguous genes at 7q11.23 commonly deleted in Williams syndrome, a multisystem developmental disorder. This gene consists of at least 14 exons, and its alternative splicing generates 3 transcript variants, all encoding the same protein. [provided by RefSeq, Jul 2008],