

## RGS19 Polyclonal Antibody

Cat No: HR1AP10086

For research use only

### Overview

Product Name	RGS19 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	23kDa
GeneID?Human?	RGS19 GAIP GNAI3IP
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>regulator of G-protein signaling 19(RGS19) Homo sapiens G proteins mediate a number of cellular processes. The protein encoded by this gene belongs to the RGS (regulators of G-protein signaling) family and specifically interacts with G protein, GAI3. This protein is a guanosine triphosphatase-activating protein that functions to down-regulate Galpha i/Galpha q-linked signaling. Alternatively spliced transcript variants encoding the same protein isoform have been found for this gene. [provided by RefSeq, Jul 2008],</p>