

KIZ Polyclonal Antibody

Cat No: HR1AP10087

For research use only

Overview

Product Name	KIZ Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	74kDa
GeneID?Human?	PLK1S1 C20orf19 KIZ NCRNA00153 HT013
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	kizuna centrosomal protein(KIZ) Homo sapiens The protein encoded by this gene localizes to centrosomes, strengthening and stabilizing the pericentriolar region prior to spindle formation. The encoded protein usually remains with the mother centrosome after centrosomal duplication. Sevral transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Feb 2013],