

DPYL1 Polyclonal Antibody

Cat No: HR1AP10083

For research use only

Overview

Product Name	DPYL1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	62kDa
GeneID?Human?	CRMP1 DPYSL1 ULIP3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>collapsin response mediator protein 1(CRMP1) Homo sapiens This gene encodes a member of a family of cytosolic phosphoproteins expressed exclusively in the nervous system. The encoded protein is thought to be a part of the semaphorin signal transduction pathway implicated in semaphorin-induced growth cone collapse during neural development. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2008],</p>