

CLSPN Polyclonal Antibody

Cat No: HR1AP10082

For research use only

Overview

Product Name	CLSPN Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	147kDa
GeneID?Human?	CLSPN
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>claspin(CLSPN) Homo sapiens The product of this gene is an essential upstream regulator of checkpoint kinase 1 and triggers a checkpoint arrest of the cell cycle in response to replicative stress or DNA damage. The protein is also required for efficient DNA replication during a normal S phase. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jun 2010],</p>