

BORA Polyclonal Antibody

Cat No: HR1AP10081

For research use only

Overview

Product Name	BORA Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	61kDa
GeneID?Human?	BORA C13orf34
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	bora, aurora kinase A activator(BORA) Homo sapiens BORA is an activator of the protein kinase Aurora A (AURKA; MIM 603072), which is required for centrosome maturation, spindle assembly, and asymmetric protein localization during mitosis (Hutterer et al., 2006 [PubMed 16890155]).[supplied by OMIM, Mar 2008],