

DAG1 Polyclonal Antibody

Cat No: HR1AP10076

For research use only

Overview

Product Name	DAG1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	98kDa
GeneID?Human?	DAG1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>dystroglycan 1(DAG1) Homo sapiens This gene encodes dystroglycan, a central component of dystrophin-glycoprotein complex that links the extracellular matrix and the cytoskeleton in the skeletal muscle. The encoded preproprotein undergoes O- and N-glycosylation, and proteolytic processing to generate alpha and beta subunits. Certain mutations in this gene are known to cause distinct forms of muscular dystrophy. Alternative splicing results in multiple transcript variants, all encoding the same protein. [provided by RefSeq, Nov 2015].</p>