

# EIF2BE Polyclonal Antibody

Cat No: HR1AP10071

For research use only

## Overview

Product Name	EIF2BE Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	79kDa
GeneID?Human?	EIF2B5 EIF2BE
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	eukaryotic translation initiation factor 2B subunit epsilon(EIF2B5) Homo sapiens This gene encodes one of five subunits of eukaryotic translation initiation factor 2B (EIF2B), a GTP exchange factor for eukaryotic initiation factor 2 and an essential regulator for protein synthesis. Mutations in this gene and the genes encoding other EIF2B subunits have been associated with leukoencephalopathy with vanishing white matter. [provided by RefSeq, Nov 2009],