

# ROS1 Polyclonal Antibody

Cat No: HR1AP10070

For research use only

## Overview

Product Name	ROS1 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	258kDa
GeneID?Human?	ROS1 MCF3 ROS
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>ROS proto-oncogene 1, receptor tyrosine kinase(ROS1) Homo sapiens This proto-oncogene, highly-expressed in a variety of tumor cell lines, belongs to the sevenless subfamily of tyrosine kinase insulin receptor genes. The protein encoded by this gene is a type I integral membrane protein with tyrosine kinase activity. The protein may function as a growth or differentiation factor receptor. [provided by RefSeq, Jul 2008],</p>