

PASK Polyclonal Antibody

Cat No: HR1AP10068

For research use only

Overview

Product Name	PASK Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	145kDa
GeneID?Human?	PASK KIAA0135
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>PAS domain containing serine/threonine kinase(PASK) Homo sapiens This gene encodes a member of the serine/threonine kinase family that contains two PAS domains. Expression of this gene is regulated by glucose, and the encoded protein plays a role in the regulation of insulin gene expression. Downregulation of this gene may play a role in type 2 diabetes. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Nov 2011],</p>