

DYRK2 Polyclonal Antibody

Cat No: HR1AP10063

For research use only

Overview

Product Name	DYRK2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	66kDa
GeneID?Human?	DYRK2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>dual specificity tyrosine phosphorylation regulated kinase 2(DYRK2) Homo sapiens DYRK2 belongs to a family of protein kinases whose members are presumed to be involved in cellular growth and/or development. The family is defined by structural similarity of their kinase domains and their capability to autophosphorylate on tyrosine residues. DYRK2 has demonstrated tyrosine autophosphorylation and catalyzed phosphorylation of histones H3 and H2B in vitro. Two isoforms of DYRK2 have been isolated. The predominant isoform, isoform 1, lacks a 5' terminal insert. [provided by RefSeq, Jul 2008],</p>