

# AHNK Polyclonal Antibody

Cat No: HR1AP10058

For research use only

## Overview

Product Name	AHNK Polyclonal Antibody
Source	Rabbit
Applications	IHC-p
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	647kDa
GeneID?Human?	AHNAK PM227
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:May be required for neuronal cell differentiation.,similarity:Contains 1 PDZ (DHR) domain.,subunit:Interacts with DYSF; the interaction is direct and Ca(2+)-independent.,