

Tmem72 Polyclonal Antibody

Cat No: HR1AP10056

For research use only

Overview

Product Name	Tmem72 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	30kDa
GeneID?Human?	TMEM72 C10orf127 KSP37
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	Transmembrane protein 72 (Kidney-specific secretory protein of 37 kDa)
Background	transmembrane protein 72(TMEM72) Homo sapiens This gene encodes a transmembrane protein which may be expressed specifically in the kidney. [provided by RefSeq, Sep 2016],