

AFG3L2 Polyclonal Antibody

Cat No: HR1AP9977

For research use only

Overview

Product Name	AFG3L2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	88kDa
GeneID?Human?	AFG3L2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	AFG3-like protein 2 (EC 3.4.24.-) (Paraplegin-like protein)
Background	AFG3 like matrix AAA peptidase subunit 2(AFG3L2) Homo sapiens This gene encodes a protein localized in mitochondria and closely related to paraplegin. The paraplegin gene is responsible for an autosomal recessive form of hereditary spastic paraplegia. This gene is a candidate gene for other hereditary spastic paraplegias or neurodegenerative disorders. [provided by RefSeq, Jul 2008],