

PON1 Polyclonal Antibody

Cat No: HR1AP9852

For research use only

Overview

Product Name	PON1 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	PON1 PON
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	Serum paraoxonase/arylesterase 1 (PON 1) (EC 3.1.1.2) (EC 3.1.1.81) (EC 3.1.8.1) (Aromatic esterase 1) (A-esterase 1) (K-45) (Serum aryldialkylphosphatase 1)
Background	paraoxonase 1(PON1) Homo sapiens The enzyme encoded by this gene is an arylesterase that mainly hydrolyzes paraoxon to produce p-nitrophenol. Paraoxon is an organophosphorus anticholinesterase compound that is produced in vivo by oxidation of the insecticide parathion. Polymorphisms in this gene are a risk factor in coronary artery disease. The gene is found in a cluster of three related paraoxonase genes at 7q21.3. [provided by RefSeq, Oct 2008].