

IgD Polyclonal Antibody

Cat No: HR1AP9848

For research use only

Overview

Product Name	IgD Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	IGHD
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	Ig delta chain C region
Background	similarity:Contains 3 Ig-like (immunoglobulin-like) domains.,