

CHID1 Polyclonal Antibody

Cat No: HR1AP9660

For research use only

Overview

Product Name	CHID1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	45kDa
GeneID?Human?	CHID1 GL008 PSEC0104 SB139
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CHID1 GL008 PSEC0104 SB139
Background	induction:Up-regulated by interleukin-4 and dexamethasone in the macrophages.,similarity:Belongs to the glycosyl hydrolase 18 family.,subunit:Interacts with STAB1.,tissue specificity:Expressed in cells of monocytic, T, B and epithelial origin.,