

Collagen alpha-1(XXVIII)Polyclonal Antibody

Cat No: HR1AP9613

For research use only

Overview

Product Name	Collagen alpha-1(XXVIII)Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	116kDa
GeneID?Human?	COL28A1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	Collagen alpha-1(XXVIII) chain
Background	collagen type XXVIII alpha 1 chain(COL28A1) Homo sapiens COL28A1 belongs to a class of collagens containing von Willebrand factor (VWF; MIM 613160) type A (VWFA) domains (Veit et al., 2006 [PubMed 16330543]).[supplied by OMIM, Nov 2010],