

AMBRA1 Polyclonal Antibody

Cat No: HR1AP9591

For research use only

Overview

Product Name	AMBRA1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	140kDa
GeneID?Human?	AMBRA1 KIAA1736
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	autophagy/beclin-1 regulator 1
Background	function:Regulates autophagy and development of the nervous system. Involved in autophagy in controlling protein turnover during neuronal development, and in regulating normal cell survival and proliferation.,similarity:Contains 3 WD repeats.,subunit:Interacts with BECN1. Probably forms a complex with BECN1 and PIK3C3.,