

Atg9a Polyclonal Antibody

Cat No: HR1AP9587

For research use only

Overview

Product Name	Atg9a Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	85kDa
GeneID?Human?	ATG9A APOG9L1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ATG9 autophagy related 9 homolog A (S. cerevisiae)
Background	function:Plays a role in autophagy.,similarity:Belongs to the ATG9 family.,