

## PGLYRP1 Polyclonal Antibody

Cat No: HR1AP9520

For research use only

### Overview

Product Name	PGLYRP1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	21kDa
GeneID?Human?	PGLYRP1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	PGLYRP1; PGLYRP; PGRP; TNFSF3L; SBBI68; Peptidoglycan recognition protein 1; Peptidoglycan recognition protein short; PGRP-S
Background	PGLYRP1 (Peptidoglycan Recognition Protein 1) is a Protein Coding gene. Among its related pathways are TNF Signaling (sino). GO annotations related to this gene include peptidoglycan binding and N-acetylmuramoyl-L-alanine amidase activity. An important paralog of this gene is PGLYRP4.