

## IL-17F Polyclonal Antibody

Cat No: HR1AP9916

For research use only

### Overview

Product Name	IL-17F Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	IL17F IL24
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	Interleukin-17F (IL-17F) (Cytokine ML-1) (Interleukin-24) (IL-24)
Background	interleukin 17F(IL17F) Homo sapiens The protein encoded by this gene is a cytokine that shares sequence similarity with IL17. This cytokine is expressed by activated T cells, and has been shown to stimulate the production of several other cytokines, including IL6, IL8, and CSF2/GM-CSF. This cytokine is also found to inhibit the angiogenesis of endothelial cells and induce endothelial cells to produce IL2, TGFB1/TGFB, and monocyte chemoattractant protein-1. [provided by RefSeq, Jul 2008],