

CD52 Polyclonal Antibody

Cat No: HR1AP9898

For research use only

Overview

Product Name	CD52 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	CD52 CDW52 HE5
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CAMPATH-1 antigen (CDw52) (Cambridge pathology 1 antigen) (Epididymal secretory protein E5) (Human epididymis-specific protein 5) (He5) (CD antigen CD52)
Background	function:May play a role in carrying and orienting carbohydrate, as well as having a more specific role.,