

DLEU7 Polyclonal Antibody

Cat No: HR1AP9855

For research use only

Overview

Product Name	DLEU7 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	DLEU7 LEU7
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	Leukemia-associated protein 7 (Deleted in lymphocytic leukemia 7)
Background	alternative products:Additional isoforms seem to exist,developmental stage:Expressed in the fetal brain, liver and kidney.,tissue specificity:Ubiquitous.,