

GSC Polyclonal Antibody

Cat No: HR1AP9863

For research use only

Overview

Product Name	GSC Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	GSC
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	Homeobox protein gooseoid
Background	gooseoid homeobox(GSC) Homo sapiens This gene encodes a member of the bicoid subfamily of the paired (PRD) homeobox family of proteins. The encoded protein acts as a transcription factor and may be autoregulatory. A similar protein in mice plays a role in craniofacial and rib cage development during embryogenesis. [provided by RefSeq, Jul 2008],