

## GMF-? Polyclonal Antibody

Cat No: HR1AP9836

For research use only

### Overview

Product Name	GMF-? Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	GMFG
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	Glia maturation factor gamma (GMF-gamma)
Background	similarity:Belongs to the actin-binding proteins ADF family. GMF subfamily.,similarity:Contains 1 ADF-H domain.,tissue specificity:Expressed predominantly in lung, heart and placenta.,