

CD238 Polyclonal Antibody

Cat No: HR1AP9577

For research use only

Overview

Product Name	CD238 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	82kDa
GeneID?Human?	KEL
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	Kell blood group, metallo-endopeptidase
Background	Kell blood group, metallo-endopeptidase(KEL) Homo sapiens This gene encodes a type II transmembrane glycoprotein that is the highly polymorphic Kell blood group antigen. The Kell glycoprotein links via a single disulfide bond to the XK membrane protein that carries the Kx antigen. The encoded protein contains sequence and structural similarity to members of the neprilysin (M13) family of zinc endopeptidases. [provided by RefSeq, Jul 2008],