

MYH6/MYH7 Polyclonal Antibody

Cat No: HR1AP9922

For research use only

Overview

Product Name	MYH6/MYH7 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	MYH6/7 MYHCA/B
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>myosin heavy chain 6(MYH6) Homo sapiens Cardiac muscle myosin is a hexamer consisting of two heavy chain subunits, two light chain subunits, and two regulatory subunits. This gene encodes the alpha heavy chain subunit of cardiac myosin. The gene is located ~4kb downstream of the gene encoding the beta heavy chain subunit of cardiac myosin. Mutations in this gene cause familial hypertrophic cardiomyopathy and atrial septal defect 3. [provided by RefSeq, Mar 2010],</p>