

## IL-17D Polyclonal Antibody

Cat No: HR1AP9917

For research use only

### Overview

Product Name	IL-17D Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	IL17D IL27 UNQ3096/PRO21175
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	Interleukin-17D (IL-17D) (Interleukin-27) (IL-27)
Background	interleukin 17D(IL17D) Homo sapiens The protein encoded by this gene is a cytokine that shares the sequence similarity with IL17. The treatment of endothelial cells with this cytokine has been shown to stimulate the production of other cytokines including IL6, IL8 and CSF2/ GM-CSF. The increased expression of IL8 induced by this cytokine was found to be NF-kappa B-dependent. [provided by RefSeq, Jul 2008],