

HSP 75 Polyclonal Antibody

Cat No: HR1AP9859

For research use only

Overview

| | |
|-----------------------|---|
| Product Name | HSP 75 Polyclonal Antibody |
| Source | Rabbit |
| Applications | IHC-p,ELISA |
| Species Reactivity | Human |
| Recommended Dilutions | |
| Immunogen | |
| Species | Rabbit |
| Storage | -20°C/1 year |
| Isotype | |
| Clonality | |
| Concentration | 1 mg/ml |
| Observed band | kDa |
| GeneID?Human? | TRAP1 HSP75 |
| Human Swiss-Prot No. | |
| Cellular localization | |
| Alternative Names | Heat shock protein 75 kDa, mitochondrial (HSP 75) (TNFR-associated protein 1) (Tumor necrosis factor type 1 receptor-associated protein) (TRAP-1) |
| Background | TNF receptor associated protein 1(TRAP1) Homo sapiens This gene encodes a mitochondrial chaperone protein that is member of the heat shock protein 90 (HSP90) family. The encoded protein has ATPase activity and interacts with tumor necrosis factor type I. This protein may function in regulating cellular stress responses. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Jan 2013], |