

# CNPY3 Polyclonal Antibody

Cat No: HR1AP9920

For research use only

## Overview

Product Name	CNPY3 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	CNPY3 CTG4A ERDA5 PRAT4A TNRC5 HSPC084 UNQ1934/PRO4409
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	Protein canopy homolog 3 (CTG repeat protein 4a) (Expanded repeat-domain protein CAG/CTG 5) (Protein associated with TLR4) (Trinucleotide repeat-containing gene 5 protein)
Background	canopy FGF signaling regulator 3(CNPY3) Homo sapiens This gene encodes a protein that binds members of the toll-like receptor protein family and functions as a chaperone to aid in folding and export of these proteins. Alternative splicing results in multiple transcript variants. Naturally occurring readthrough transcription occurs between this locus and the downstream GNMT (glycine N-methyltransferase) gene and is represented with GeneID:107080644. [provided by RefSeq, Jan 2016],