

EYA1/EYA4 Polyclonal Antibody

Cat No: HR1AP9919

For research use only

Overview

Product Name	EYA1/EYA4 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	EYA3/4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>EYA transcriptional coactivator and phosphatase 1(EYA1) Homo sapiens This gene encodes a member of the eyes absent (EYA) family of proteins. The encoded protein may play a role in the developing kidney, branchial arches, eye, and ear. Mutations of this gene have been associated with branchiootorenal dysplasia syndrome, branchiootic syndrome, and sporadic cases of congenital cataracts and ocular anterior segment anomalies. A similar protein in mice can act as a transcriptional activator. Alternatively spliced transcript variants have been identified for this gene. [provided by RefSeq, Dec 2013],</p>