

ERLIN1/2 Polyclonal Antibody

Cat No: HR1AP9618

For research use only

Overview

Product Name	ERLIN1/2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	40+42kDa
GeneID?Human?	ERLIN1/ERLIN2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	similarity:Belongs to the band 7/mec-2 family.,subcellular location:Associated with lipid raft-like domains of the endoplasmic reticulum membrane.,tissue specificity:Expressed in heart, placenta, liver, kidney, pancreas, prostate, testis, ovary and small intestin.,