

Imp3 Polyclonal Antibody

Cat No: HR1AP9642

For research use only

Overview

Product Name	Imp3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	25kDa
GeneID?Human?	IMP3 C15orf12 MRPS4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	IMP3 C15orf12 MRPS4
Background	IMP3, U3 small nucleolar ribonucleoprotein(IMP3) Homo sapiens This gene encodes the human homolog of the yeast Imp3 protein. The protein localizes to the nucleoli and interacts with the U3 snoRNP complex. The protein contains an S4 domain. [provided by RefSeq, Jul 2008],