

Nanos Homologue 1 (NANOS1) Polyclonal Antibody

Cat No: HR1AP9617

For research use only

Overview

Product Name	Nanos Homologue 1 (NANOS1) Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	30kDa
GeneID?Human?	NANOS1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	Nanos homolog 1 (NOS-1) (EC_Rep1a)
Background	nanos C2HC-type zinc finger 1(NANOS1) Homo sapiens This gene encodes a CCHC-type zinc finger protein that is a member of the nanos family. This protein co-localizes with the RNA-binding protein pumilio RNA-binding family member 2 and may be involved in regulating translation as a post-transcriptional repressor. Mutations in this gene are associated with spermatogenic impairment. [provided by RefSeq, Sep 2015],