

SCRIB Polyclonal Antibody

Cat No: HR1AP12998

For research use only

Overview

Product Name	SCRIB Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	179kDa
GeneID?Human?	SCRIB CRIB1 KIAA0147 LAP4 SCRB1 VARTUL
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>scribbled planar cell polarity protein(SCRIB) Homo sapiens This gene encodes a protein that was identified as being similar to the Drosophila scribble protein. The mammalian protein is involved in tumor suppression pathways. As a scaffold protein involved in cell polarization processes, this protein binds to many other proteins. The encoded protein binds to papillomavirus E6 protein via its PDZ domain and the C-terminus of E6. Two alternatively spliced transcript variants that encode different protein isoforms have been found for this gene. [provided by RefSeq, Nov 2011],</p>