

TWST1 Polyclonal Antibody

Cat No: HR1AP12980

For research use only

Overview

Product Name	TWST1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	22kDa
GeneID?Human?	TWIST1 BHLHA38 TWIST
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>twist family bHLH transcription factor 1(TWIST1) Homo sapiens Basic helix-loop-helix (bHLH) transcription factors have been implicated in cell lineage determination and differentiation. The protein encoded by this gene is a bHLH transcription factor and shares similarity with another bHLH transcription factor, Dermo1. The strongest expression of this mRNA is in placental tissue; in adults, mesodermally derived tissues express this mRNA preferentially. Mutations in this gene have been found in patients with Saethre-Chotzen syndrome. [provided by RefSeq, Jul 2008].</p>