

## DRAM1 Polyclonal Antibody

Cat No: HR1AP12547

For research use only

### Overview

Product Name	DRAM1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	26kDa
GeneID?Human?	DRAM1 DRAM
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	DNA damage regulated autophagy modulator 1(DRAM1) Homo sapiens This gene is regulated as part of the p53 tumor suppressor pathway. The gene encodes a lysosomal membrane protein that is required for the induction of autophagy by the pathway. Decreased transcriptional expression of this gene is associated with various tumors. This gene has a pseudogene on chromosome 4. [provided by RefSeq, Jul 2008],