

GPR84 Polyclonal Antibody

Cat No: HR1AP12594

For research use only

Overview

Product Name	GPR84 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	43kDa
GeneID?Human?	GPR84 EX33
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:Orphan receptor.,similarity:Belongs to the G-protein coupled receptor 1 family.,tissue specificity:Expressed in brain, heart, muscle, colon, thymus, spleen, kidney, liver, placenta, intestine, bone marrow, lung and peripheral blood leukocytes. Within leukocyte population more abundant expression in neutrophils and eosinophils, in comparison with that in T- or B-lymphocytes.,