

DLL3 Polyclonal Antibody

Cat No: HR1AP12418

For research use only

Overview

Product Name	DLL3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	67kDa
GeneID?Human?	DLL3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>delta like canonical Notch ligand 3(DLL3) Homo sapiens This gene encodes a member of the delta protein ligand family. This family functions as Notch ligands that are characterized by a DSL domain, EGF repeats, and a transmembrane domain. Mutations in this gene cause autosomal recessive spondylocostal dysostosis 1. Two transcript variants encoding distinct isoforms have been identified for this gene. [provided by RefSeq, Jul 2008],</p>