

## NR1H4 Polyclonal Antibody

Cat No: HR1AP12218

For research use only

### Overview

Product Name	NR1H4 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	53kDa
GeneID?Human?	NR1H4 BAR FXR HRR1 RIP14
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	nuclear receptor subfamily 1 group H member 4(NR1H4) Homo sapiens This gene encodes a ligand-activated transcription factor that shares structural features in common with nuclear hormone receptor family members. This protein functions as a receptor for bile acids, and when bound to bile acids, binds to DNA and regulates the expression of genes involved in bile acid synthesis and transport. Alternatively spliced transcript variants encoding different isoforms have been described. [provided by RefSeq, Feb 2016],