

JAG2 Polyclonal Antibody

Cat No: HR1AP12151

For research use only

Overview

Product Name	JAG2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	136kDa
GeneID?Human?	JAG2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>jagged 2(JAG2) Homo sapiens The Notch signaling pathway is an intercellular signaling mechanism that is essential for proper embryonic development. Members of the Notch gene family encode transmembrane receptors that are critical for various cell fate decisions. The protein encoded by this gene is one of several ligands that activate Notch and related receptors. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008].</p>