

NUFP2 Polyclonal Antibody

Cat No: HR1AP11706

For research use only

Overview

Product Name	NUFP2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	76kDa
GeneID?Human?	NUFIP2 KIAA1321 PIG1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function: Binds RNA.,PTM: Phosphorylated upon DNA damage, probably by ATM or ATR.,subcellular location: Distribution is cell cycle-modulated, being cytoplasmic in the G2/M phase and accumulating in nucleus during the G1 phase.,subunit: Interacts with FMRP.,