

M3K14 Polyclonal Antibody

Cat No: HR1AP11651

For research use only

Overview

Product Name	M3K14 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	104kDa
GeneID?Human?	MAP3K14 NIK
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>mitogen-activated protein kinase kinase kinase 14(MAP3K14) Homo sapiens This gene encodes mitogen-activated protein kinase kinase kinase 14, which is a serine/threonine protein-kinase. This kinase binds to TRAF2 and stimulates NF-kappaB activity. It shares sequence similarity with several other MAPKK kinases. It participates in an NF-kappaB-inducing signalling cascade common to receptors of the tumour-necrosis/nerve-growth factor (TNF/NGF) family and to the interleukin-1 type-I receptor. [provided by RefSeq, Jul 2008].</p>