

## CENPN Polyclonal Antibody

Cat No: HR1AP11637

For research use only

### Overview

Product Name	CENPN Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	37kDa
GeneID?Human?	CENPN C16orf60 ICEN32 BM-309
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	centromere protein N(CENPN) Homo sapiens The protein encoded by this gene forms part of the nucleosome-associated complex and is important for kinetochore assembly. It is bound to kinetochores during S phase and G2 and recruits other proteins to the centromere. Pseudogenes of this gene are located on chromosome 2. Alternative splicing results in multiple transcript variants that encode different protein isoforms. [provided by RefSeq, Jul 2012].