

## ZO1 Polyclonal Antibody

Cat No: HR1AP11467

For research use only

### Overview

Product Name	ZO1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	192kDa
GeneID?Human?	TJP1 ZO1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	tight junction protein 1(TJP1) Homo sapiens This gene encodes a protein located on a cytoplasmic membrane surface of intercellular tight junctions. The encoded protein may be involved in signal transduction at cell-cell junctions. Alternative splicing of this gene results in multiple transcript variants. [provided by RefSeq, Jul 2014],