



SFTPD Polyclonal Antibody

Cat No: HR1AP11183

For research use only

Overview

Product Name	SFTPD Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	41kDa
GeneID?Human?	SFTPD COLEC7 PSPD SFTP4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	surfactant protein D(SFTPD) Homo sapiens The protein encoded by this gene is part of the innate immune response, protecting the lungs against inhaled microorganisms and chemicals. The encoded protein may also be involved in surfactant metabolism. [provided by RefSeq, Jul 2015],