

## GNAO Polyclonal Antibody

Cat No: HR1AP10813

For research use only

### Overview

Product Name	GNAO Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	38kDa
GeneID?Human?	GNAO1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	G protein subunit alpha o1(GNAO1) Homo sapiens The protein encoded by this gene represents the alpha subunit of the Go heterotrimeric G-protein signal-transducing complex. Defects in this gene are a cause of early-onset epileptic encephalopathy. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Aug 2015],